

Sending Digital Images to RCS

Follow the steps below to send your digital images to the RCS faster, safer and cheaper.

Step One

Make sure that your digital equipment can produce DICOM or DICOM-compatible images.

Step Two

Create a secure tunnel that allows you to send confidential information over an internet connection which is fully HIPAA compliant. Your IT support staff will work with Dan Prill in our IT department to create this secure tunnel. You will need a firewall or router that supports IPSEC tunnels.

Step Three

Contact your PACS vendor who will specify certain parameters in your DICOM software to be able to send the images to our server.

These settings are as follows:

AE Title: MN_NC_AC_STORE
IP Address: 172.16.1.29
Port: 4321

Step Four

Fax an RCS transmittal form to 952-886-7565 each time you send a case. The fax indicates to us that your films are on our work list to be read.

Sending images on CDs

If your digital images are not DICOM or DICOM-compatible you can send them in on a CD. Get a transmittal form from our website, fill it out and send it in with your CD. We typically don't return CDs unless they are third-party generated or you request that they be returned.

For assistance call

Anita L. Manne, DC, DACBR, Director of RCS
952-885-5464

Dan Prill, Information Technology
952-887-1399